

HIGHER PERFORMANCE FOR HIGHER PERFORMING CHAIRS

Introducing the Interstate Batteries® DCL series



DCL-55



DCL-75



HIGH PERFORMANCE FOR HIGH-END EQUIPMENT

Today's top-tier power wheelchairs need a higher-performing battery. Interstate® has your customers covered with our new DCL series. These new batteries are designed for better charge acceptance with lower-voltage chargers. A better charge leads to happier customers.




The Interstate® DCL-55 and DCL-75 batteries feature:

- Better Performance With Low-voltage Chargers
- Deep Cycle Performance With a Nonspillable Design
- Sealed, Absorbed Glass-mat Technology
- Travel Certification
- No Break-in Cycle
- 1-year Free-replacement Warranty

| PART NUMBER | INDUSTRY SIZE | RATING | HCPCS* CODE | DIMENSIONS L x W x H | TERMINAL AVAILABILITY | TERMINAL TORQUE |
|-------------|---------------|--------|-------------|----------------------|--|-----------------|
| DCL-55 | 22NF | 55 Ah | E2361 | 8.98 x 5.39 x 9.06 | 1/4" -20 Insert or U (Hardware Included) | 30-50 in-lbs. |
| DCL-75 | 24 | 75 Ah | E2363 | 10.2 x 6.61 x 8.43 | 1/4" -20 Insert or U (Hardware Included) | 100-130 in-lbs. |

Charging Voltage: Float Charge Voltage 13.5 to 13.8 VDC at 77° F (25° C), Cycle Charge Voltage 14.3 to 15.0 VDC at 77° F (25° C).



U TERMINAL



* HCPCS codes are subject to change and are shown here for power chair applications only. Before submitting billing to Medicare, please contact the DMERC for verification of codes. Interstate Batteries disclaims all liability and responsibility for improper use of HCPCS codes by consumers, resellers or other end users.